

Real-time plotter and fitter server manual

Oleksiy Onishchenko

March 20, 2021

1 Notation

- File and folder names are given in the following font *myfavoritefile.py*
- Angled brackets with a datatype inside, so like `<integer>` mean that the user must provide one instance of the corresponding data type. It is possible that one has to provide lists, in which case the data type of the entries in the list will also be provided in angles brackets, so say a list of two strings would be `<list(<string>,<string>)>`. If there are several data types possible (which should not happen too much), they will be separated by a comma (so say `<integer, string>`); do NOT interpret this as a list of two entries, this is integer OR string, one set of single brackets means only one instance. If we don't want to specify the datatype, we will use the word "any", and we will use ... to denote an undefined number; for example, a list of integers of arbitrary length will be given as `<list(<integer>,...)>`
- Three

2 Basic principles, design philosophy, and governing ideas

3 Client-server structure

The principal idea of the project is to make this plotter-fitter work as a remote server. It should be possible to send all data and commands remotely and get answers back. Of course there is a GUI layer on top, in order to visualize the plots and fits, read the results off the screen, and also do some manual adjustment of plotting, but a big focus of programming this tool is to make it function as a remote server.

The main server class `TCPIPserver` is located in file *socketserver.py*. It is defined to work in two modes: either one-way, where it only receives messages from the client and does not send any info back, or in two-way more, where a round of communication consists of message reception and transmission. ~~At some point in the future this should probably be combined into a single server without a real distinction in the functions themselves whether~~

4 How the fitter itself functions

TCIPserver calls class `GeneralFitter1D` with an instance of `Fitmodel` class as the only parameter. Fitting itself is done in function `GeneralFitter1D.doFit()`, meaning that the optimizer from (scipy) is called in that function.

5 TCP/IP commands

5.1 General format of commands

All commands must be sent as character strings in **utf-8** encoding (if that's impossible, one could implement ascii encoding/decoding procedure, but it's not done yet). Each individual command must conform to JSON-RPC2.0 standard. Therefore the strings look as follows:

JSON-RPC2.0 command format

```
{ "jsonrpc": "2.0", "method": <string>, "params":  
<dictionary>, "id": <integer> }
```

where the dictionary corresponding to "params" is a data structure corresponding to the python dictionary and has the format {<string>:<any>,...}, so it is a list of comma-separated pairs, of arbitrary length, where inside each pair itself the entries are separated by a colon. The first element of each pair is a string (that's known as the *key*), and the second element can be in principle any datatype, including a dictionary itself (that's known as *value*). Note the curly braces around: they must be there.

For those who program in Python and understand the lingo: the value corresponding "params" has the standard form of a Python dictionary, with all keys being strings. For those who do not program in Python and do not understand the lingo: make sure to build this dictionary exactly in the format described above.

5.2 Available "method" values and responses from the server

As of now, the implemented methods are

- `'doClear'` This will clear the information from the screen and from the server memory.
- `'setConfig'` This will send configuration parameters to the plotter and fitter, such as axis labels, plot label, legend labels, etc
- `'addData'` This send data point by point to the plotter, or alternatively (later) lists of data points at once. Once the plotter receives each data point, it immediately puts it on the screen

- `'doFit'` This will perform the fit to (some of) the data that has been previously sent to the server. Cropping is also done here, because cropping is something that's used only for fitting.
- `'getFitResult'` This tells the plotter which fit to send back to the client.
- `'getConfig'` Not implemented yet, but envisioned to get configurations back to the server

whenever the method is anything other than “getFitResult” or “getConfig”, the server send to the client a JSON-RPC2.0 string of the following form:

```
{ 'jsonrpc': '2.0', 'result': 'MessageReceived' }
```

Whenever the method is “getFitResult”, the server will respond as follows:

```
{ 'jsonrpc': '2.0', 'result': <dictionary> }
```

and the result dictionary will contain fit results in the form like

```
{ 'frequency': { 'fitvalue': 100000., 'fitererror': 3000. },
... , 'costfunction': 200. }
```

so the result will contain a dictionary of fit parameters with their corresponding fit values and possibly fit errors, and finally it will also contain the minimum cost function that was obtained after optimization. One could also imagine extending this dictionary of results to give more information about the fit, but that should be easy to do given that this is simply extending the dictionary, and one has the tree structure such that on the right inside each colon-separated pair one can have another dictionary (Python lingo: dictionaries can have dictionaries as their values).

5.3 Available “params” values

As we have seen, “params” in the request must be sent as dictionaries, so in the form {<string>:<any>,...}. The following table summarizes what can go into these dictionaries

Sending “params” to the server

We list the parameters with an explanation of the possible values (the parameter literal is before the colon, the possible values are after the colon, explanation in parentheses

Case 1: “method” is “doClear”

- “clearData” : “all” (This clears everything, so data, fits, axis labels, etc.)
- “clearData”: <int> (This clears the curve, data and corresponding fit if available, specified by the integer)
- “clearConfig”: “all” (This clears the axis labels and the title)

For later implementation: clearing individual curves, clearing individual labels, etc.

Case 2: “method” is “setConfig”

- “setAxisLabels” : <list(<string>,<string>)> (first one in the list is x-axis label, second one is y-axis label)
- “setPlotTitle” : <string>
- “setPlotLegend” : <dictionary(<string>:<string>,...)> (The first string in the dictionary has form “curve1” for example, etc. which just defines the curve to use, the second string will be the label that we want to put in the legend for that curve)

... (more commands can be added later)

Case 3: “method” is “addData”

- Adding a single data point (so a point to a single curve)

“dataPoint” : <dictionary>.

The dictionary that goes with the dataPoint have the form

```
{ 'curveNumber':<integer>, 'xval':<float,int>,  
'yval':<float,int>, 'yerr':<float,int>,  
'xerr':<float,int> },
```

and in this case the options “xerr” and “yerr” are not required, they are optional. This will add a single data point to the curve that has been specified.

- Adding data points to multiple curves at the same time

“pointList” : <list(<dictionary>,...)>.

and each dictionary has exactly the same format as in the item above. The list is specified by using these symbols before and after:

```
[ ]
```

so it's formatted as a standard Python list.

This will simply multiple points (to different curves) in the same round of JSON-RPC2.0 communication.

Notice that this list is parsed directly in the JSONInterpreter.py, so the main program gets directly a list of commands corresponding to “dataPoint”.

(this can be modified for adding points to curves, but the logic of it being a dictionary will stay the same)

Case 4: “method” is “doFit”

- “fitFunction”: <str> (Name of the fitting function to be used in case “performFitting” is called afterwards. This name must be defined in fitmodels.py)

- “curveNumber”: <int> (The curve that will be processed in this particular call iteration of doFit routines)

- “startingParameters” : <dictionary>

This dictionary must come in the form

{ 'frequency':<float>, 'center':<float>,...}, so these will be the names of the variables are they are defined in the fit model. This depends of course on what one wants to fit, Gaussian, sinusoidal, Lorentzian, etc.

Whatever parameters from the fit model are not specified will be handled by the automatic parameter estimation routine (which, OK, could be good or bad). This parameter is in principle optional, so one can not specify it at all, in which case all starting parameters will be handled by the automatic initial parameter estimation routine.

- “startParamsLimits” : <dictionary> (This is optional, this will tell the fitter the limits of paramter search. For simple fits, one can leave this to fitting algorithm defaults)

{ 'frequency':<list>(<float>,<float>), 'center':list(<float>,<float>),...}, one must of course use the same names for fit parameters, but then after the colon one must provide the lower and upper limit for each fit parameter. It is not required to provide the limits for each parameter, however if one does provide limits, one must provide both lower and upper limits. One cannot simply leave one limit empty.

- “cropLimits” : <list(<float, -inf>,<float,inf>)> (This will crop the data before sending it to the fitter. This item is optional, if not provided, all data for the given curve will be used. If one wants a one-sided limit, then the other side must be specified as -inf or inf)

- “curveNumber” : <integer> (This just tells which curve we are supposed to fit)

- “fitMethod” : <dictionary> (This is not yet implemented, but this will give access to the different possibilities of scipy optimizer routines)

Important note: All string on the left of the colon (so things like “clearData”, “axisLabels”, “plotTitle”, so on, must consist of a word starting with a lowercase letter, followed by one or more words starting with an uppercase letter, without spaces, so a possible dictionary to go with params could have the form like { 'firstSecond': <any> }. This is crucial. Inside the program, it does string parsing by detecting the locations of the capital letters, and then it calls the

corresponding functions, which have exactly the same name but with no capital letter and with an underscore separator. So the function called when parameter “firstSecond” will be called “first_second”. See Section 2 for a more detailed explanation of the structure.

Case 5: “method” is “getFitResult”

- “curveNumber” : <integer> (This will simply return the fit result for the chosen curve)

6 Background software and hardware requirements